


50mL Polypropylene Conical Centrifuge Tubes with Printed Graduations

Item# 6288/6289/6287



Item# 6288

Description

General purpose centrifuge tube, 500/case. 50mL capacity. Blue screwcap. Black printed graduations. White printed writing space. Available in sterile with attached caps and non-sterile with separate caps.

Features

- Tube molded from chemical resistant, hydrophobic polypropylene (PP)
- Cap molded from high density polyethylene (HDPE)
- Conical style bottom
- Printed black graduations (0.5 mL graduations)
- White writing area

Products

Item#	Description	Packaging
6288	50mL, PP, sterile, 25 per bag	500 per case
6289	50mL, PP, sterile, 25 per rack	500 per case
6287	50mL, PP, non-sterile, bulk, with separate blue screwcap	500 per case



Item# 6289 in 25 place rack

50mL Polypropylene Conical Centrifuge Tubes with Printed Graduations


Item# 6288/6289/6287

Technical Specifications

Material: Tube: Hydrophobic Polypropylene (PP)
Cap: High Density Polyethylene (HDPE)
Sterilization: Ethylene Oxide (EO)
Autoclavable: Yes (without cap)
Temperature range: -40°C to +80°C
Spin rate: 8500 RCF

Product Certifications

- ISO 13485
- CE



Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

Globe Scientific Inc.

610 Winters Avenue • Paramus, NJ 07652

www.globescientific.com

© Copyright 2016 Globe Scientific Inc. All rights reserved.