



Pipette tips Filter tips Volume range	Transferpette® single channel** Nominal volume															Transferpette® multichannel** Nominal volume							
	1µL	2.5µL	5µL	10µL	20µL	20µL***	25µL	50µL	100µL	200µL	500µL	1000µL	2mL	5mL	10mL	10µL	20µL	25µL	50µL	100µL	200µL	300µL	
	0.1 - 20µL	✓	✓		✓	✓											✓	✓					
0.5 - 20µL		✓		✓	✓											✓	✓						
1 - 50µL		✓		✓	✓											✓	✓						
2 - 200µL*			✓			✓	✓	✓	✓	✓								✓	✓	✓	✓	✓	✓
5 - 300µL			✓			✓	✓	✓	✓	✓								✓	✓	✓	✓	✓	✓
50 - 1000µL*											✓	✓											
0.5 - 5mL													✓	✓									
1 - 10mL															✓								
0.1 - 1µL	✓			✓												✓							
0.5 - 10µL		✓		✓	✓											✓	✓						
1 - 20µL		✓		✓	✓											✓	✓						
2 - 20µL			✓			✓	✓	✓	✓	✓								✓	✓	✓	✓	✓	✓
5 - 100µL			✓			✓	✓	✓	✓	✓								✓	✓	✓	✓	✓	✓
5 - 200µL										✓								✓	✓	✓	✓	✓	✓
50 - 1000µL											✓	✓											

✓ = Tip volume less than pipettes nominal volume

*) Bulk tips are yellow or blue colored, racked tips are colorless in a yellow-colored or blue-colored tip tray

***) Electronic pipettes are not available in all specified volumes *** Transferpette® with yellow color code